

Runthrough™ NS

CONTROL IS IN **YOUR HANDS**

Runthrough® NS features excellent torque and exceptional performance, even in the most challenging cases.




Runthrough™ NS


PRODUCT FEATURES

The Duo Core direct-joint technology of Runthrough NS family delivers definitive torque response and easier, more flexible steering control

- The directly jointed stainless steel and nitinol shaft provides a superior torque transmission.
- The nitinol alloy core to tip ensures a superior durable shaping memory.
- The “M Coat” hydrophilic coating enables smooth manipulation even in tortuous anatomy



PRODUCT SPECIFICATIONS



Diameter:
0.014" / 0.36 mm


Length of the radiopaque marker:
3 cm

Distal format:
Moldable
straight


Distal lining:
Hydrophilic M Coat™

Useful length:
180 cm

Runthrough™ NS has an excellent torque of 1:1 and an exceptional performance even in long-lasting procedures.



Coil of the Distal Tip



Proximal markers for a better differentiation

Runthrough™ NS Floppy:	Without mark	
Runthrough™ NS ExtraFloppy:	1 marker	
Runthrough™ NS Hypercoat:	2 markers	
Runthrough™ NS Intermediate: ...	3 markers	



PUSHING BOUNDARIES

Terumo Interventional Systems is committed to your success, providing innovative solutions and support for your most challenging cases.

We tirelessly search for new ways to help you apply **effective solutions and achieve the best results for your patients.**

Type	Product Code	Outer Diameter (inches)	Length (cm)	Load of the tip	Radio-pacity (cm)	Format of the tip	Support	Lining
Floppy	TW-AS418FA	0.014	180	1 gf	3	Moldable straight	Mild	Hydrophobic and Hydrophilic
Hypercoat	TW-AS418FH	0.014	180	1 gf	3	Moldable straight	Low	Hydrophilic
Intermediate	TW-AS418IA	0.014	180	3.6 gf	3	Moldable straight	High	Hydrophobic and Hydrophilic
Extra Floppy	TW-AS418XA	0.014	180	0.6 gf	3	Moldable straight	High	Hydrophobic and Hydrophilic

MORE INFORMATION



TERUMO AMERICAS HOLDING, INC.
2101 Cottontail Lane, Somerset,
New Jersey, 08873 U.S.A.
TEL: +1-732-302-4900
Tollfree (USA): +1-800-283-7866

TERUMO MEDICAL DE MEXICO S.A. DE C.V.
Av. Insurgentes Sur 1647,
Piso 11 Col. San José Insurgentes
Del. Benito Juárez, C.P. 03900
México, D.F., México
TEL: +52-55-1085-0770

TERUMO COLOMBIA ANDINA S.A.S.
Calle 100 #23-10 Piso 2
Bogotá D.C. Colombia Tel:
+57-1256-8400

TERUMO LATINAMERICA CORPORATION
Doral Corporate Center I,
8750 NW 36th Street, Suite 240,
Miami, FL 33178 U.S.A.
Tel: +1-305-477-4822

TERUMO MEDICAL DO BRASIL LTDA.
Praça General Gentil Falcão, No.
108-9º andar
Cidade Monções - Brooklin Novo
CEP: 04571-150
São Paulo - SP, Brasil
TEL: +55-11-3594-3800

TERUMO CHILE LTDA.
Avenida Apoquindo 2929
Piso 16 Las Condes Santiago
de Chile 7550246
TEL: +56-2-480-9600

©TERUMO CORPORATION. June, 2021.

All the names are trademarks or registered trademarks of TERUMO CORPORATION and their respective owners.

TERUMO CORPORATION
2-44-1 Hatagaya, Shibuya-Ku, Tokyo, 151-0072 Japan www.terumo.co.jp